

BLOOD PRESSURE MONITOR BP01 - GENERATION 2

Telli Health's 4G Blood Pressure Monitor is the most accurate 4G cellular device on the market. This FDA- and CE-certified device combines the traditional testing kit (with multiple components) into one simple device with cellular-enabled 4G technology. This device is easy-to-use, durable, reliable, and accurate.


KEY UPGRADES

- Improved cellular connectivity
- Improved battery management
- Sends device status to platform
- Automatically flags readings sent from memory
- Configurable reminder alert to take a reading



FEATURES

- Readings in memory automatically sent during check-ins
- Operational device data sent, i.e. signal strength and battery level
- RPM controllable Reminder Alerts
 - Start / stop time settable
 - Cleared upon reading
 - RPM can collect adherence behavior / build AI
- Reminder silencing automatically sent to RPM platform for improved patient adherence data
- One-touch operation: Just press the "Start" button once to measure your blood pressure and heart rate. The entire process takes less than 30 seconds. Backlit and largest display of all connected blood pressure monitors in the marketplace
- Voice operated - English Language
- Device memory bank that can store up to 180 readings which can be forwarded upon cellular connection
- The most accurate 4G cellular blood pressure monitor on the market with an accuracy of ± 3 mmHg.



SPECIFICATIONS

- ✓ Visual display of diastolic, systolic, and pulse rate data once the monitor has taken a reading.
- ✓ One power option: lithium-ion powered and rechargeable
- ✓ Adjustable adult-size BP cuff monitor for upper arm 22cm – 42cm / 8.7in – 16.5in circumference. The left side of the monitor screen displays a color code to distinguish blood pressure levels.
- ✓ Fully automatic voice indication
- ✓ Output power: 6V \pm 5%
- ✓ Pressure measurement: 0 – 299 mmHg (accuracy: \pm 3 mmHg)
- ✓ Pulse measurement: 40 – 199 pulses/min (accuracy: \pm 5%)
- ✓ Air pressure mode: micro air pump automatically inflates and electronic automatic control deflates
- ✓ Automatic shutdown: about 3 minutes



THIS KIT CONTAINS

- ✓ BP device
- ✓ Adjustable 22cm – 42cm cuff with tube (other cuff sizes available)
- ✓ Rechargeable lithium-ion battery
- ✓ Cable
- ✓ USB charging adapter
- ✓ Instruction manual



**tellihealth**

www.tellihealth.com